

Product Datasheet

RXRB Antibody (orb1241430)


Description

RXRB Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RXRB.
Target	RXRB
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RXRB antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	RXRB antibody can be used for detection of RXRB by ELISA at 1:62500. RXRB antibody can be used for detection of RXRB by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	57 kDa
Uniprot ID	P28702
NCBI	NP_068811
Dilution Range	RXRB antibody can be used for detection of RXRB by ELISA at 1:62500. RXRB antibody can be used for detection of RXRB by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Brain at 0.001 µg/mL.



Antibody used in WB on Hum. Adult Placenta at 0.001 µg/mL.



Antibody used in WB on Hum. Fetal Muscle at 0.001 µg/mL.