
Product Datasheet

SLC43A2 Antibody (orb1241524)

Description

SLC43A2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

 Antibody
used in WB
on Human
Jurkat 0.2-...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC43A2.

Target

SLC43A2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC43A2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SLC43A2 antibody can be used for detection of SLC43A2 by ELISA at 1:62500. SLC43A2 antibody can be used for detection of SLC43A2 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

53 kDa, 63 kDa

Uniprot ID
[Q8N370](#)
NCBI
[AAH71859](#)
Dilution Range

SLC43A2 antibody can be used for detection of SLC43A2 by ELISA at 1:62500. SLC43A2 antibody can be used for detection of SLC43A2 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.