
Product Datasheet

SLC25A25 Antibody (orb1241571)

Description

SLC25A25 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

 Antibody
used in WB
on Human
HepG2 at
0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC25A25.

Target

SLC25A25

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC25A25 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SLC25A25 antibody can be used for detection of SLC25A25 by ELISA at 1:312500. SLC25A25 antibody can be used for detection of SLC25A25 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

55 kDa, 55 kDa, 53 kDa

Uniprot ID
[Q6KCM6](#)
NCBI
[NP_001006642](#)
Dilution Range

SLC25A25 antibody can be used for detection of SLC25A25 by ELISA at 1:312500. SLC25A25 antibody can be used for detection of SLC25A25 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.