
Product Datasheet

FAM129A Antibody (orb1241680)

Description	FAM129A Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Transfected 293T ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human FAM129A.	
Target	FAM129A	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store FAM129A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	FAM129A antibody can be used for detection of FAM129A by ELISA at 1:312500. FAM129A antibody can be used for detection of FAM129A by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	103 kDa	
Uniprot ID	Q9BZQ8	
NCBI	NP_443198	
Dilution Range	FAM129A antibody can be used for detection of FAM129A by ELISA at 1:312500. FAM129A antibody can be used for detection of FAM129A by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	