

Product Datasheet

UPB1 Antibody (orb1241715)

Description	UPB1 Antibody	
Species/Host	Rabbit	
Reactivity	C. elegans, Canine, Drosophila, Human, Mouse, Plant, Rat, Zebrafish	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human Jurkat 2.5 ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human UPB1.	
Target	UPB1	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store UPB1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	UPB1 antibody can be used for detection of UPB1 by ELISA at 1:1562500. UPB1 antibody can be used for detection of UPB1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	42 kDa	
Uniprot ID	Q9UBR1	
NCBI	NP_057411	
Dilution Range	UPB1 antibody can be used for detection of UPB1 by ELISA at 1:1562500. UPB1 antibody can be used for detection of UPB1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	