

Product Datasheet

MMP1 Antibody (orb1241782)

Description

MMP1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MMP1.

Target

MMP1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store MMP1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

MMP1 antibody can be used for detection of MMP1 by ELISA at 1:312500. MMP1 antibody can be used for detection of MMP1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW


52 kDa

Uniprot ID
P03956
NCBI
NP_002412
Dilution Range


MMP1 antibody can be used for detection of MMP1 by ELISA at 1:312500. MMP1 antibody can be used for detection of MMP1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Intestine ...



Antibody used in IHC on Human Macrophage...



Antibody used in IHC on Human Intestine ...