

Product Datasheet

AOC1 Antibody (orb1241800)

Description	AOC1 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	Antibody used in WB on Human HepG2 at 1 ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ABP1.	
Target	AOC1	Antibody used in IHC on Human kidney at ...
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ABP1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	ABP1 antibody can be used for detection of ABP1 by ELISA at 1:312500. ABP1 antibody can be used for detection of ABP1 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	83 kDa	
Uniprot ID	P19801	
NCBI	NP_001082	
Dilution Range	ABP1 antibody can be used for detection of ABP1 by ELISA at 1:312500. ABP1 antibody can be used for detection of ABP1 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	