



Product Datasheet

GSTM2 Antibody (orb1241819)

Description	GSTM2 Antibody	
Species/Host	Rabbit	Antibody used in WB on Human Liver at 0...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GSTM2.	
Target	GSTM2	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GSTM2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in IHC on Human Skeletal M...
Note	For research use only	
Application notes	GSTM2 antibody can be used for detection of GSTM2 by ELISA at 1:312500. GSTM2 antibody can be used for detection of GSTM2 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	26 kDa	Antibody used in IHC on Human Skeletal m...
Uniprot ID	P28161	
NCBI	NP_000839	
Dilution Range	GSTM2 antibody can be used for detection of GSTM2 by ELISA at 1:312500. GSTM2 antibody can be used for detection of GSTM2 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	