



Product Datasheet

FTCD Antibody (orb1241866)

Description	FTCD Antibody
Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human FTCD.
Target	FTCD
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store FTCD antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	FTCD antibody can be used for detection of FTCD by ELISA at 1:62500. FTCD antibody can be used for detection of FTCD by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	60 kDa, 59 kDa
Uniprot ID	O95954
NCBI	NP_006648
Dilution Range	FTCD antibody can be used for detection of FTCD by ELISA at 1:62500. FTCD antibody can be used for detection of FTCD by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Fetal Liver...



Antibody used in IHC on Human Skin at 4....