
Product Datasheet

PRPS2 Antibody (orb1241887)

Description

PRPS2 Antibody

Species/Host

Rabbit

Reactivity

Canine, Drosophila, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested Applications

ELISA, IHC, WB

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PRPS2.

Target

PRPS2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PRPS2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PRPS2 antibody can be used for detection of PRPS2 by ELISA at 1:312500. PRPS2 antibody can be used for detection of PRPS2 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

35 kDa, 35 kDa, 35 kDa

Uniprot ID
[P11908](#)
NCBI
[NP_002756](#)
Dilution Range


PRPS2 antibody can be used for detection of PRPS2 by ELISA at 1:312500. PRPS2 antibody can be used for detection of PRPS2 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human HepG2 at 0....



Antibody used in IHC on Human Kidney at ...