

Product Datasheet

HNRNPA3 Antibody (orb1241936)

Description

HNRNPA3 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human HNRPA3.

Target

HNRNPA3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store HNRPA3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

HNRPA3 antibody can be used for detection of HNRPA3 by ELISA at 1:62500. HNRPA3 antibody can be used for detection of HNRPA3 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW


42 kDa

Uniprot ID
[P51991](#)
NCBI
[NP_919223](#)
Dilution Range


HNRPA3 antibody can be used for detection of HNRPA3 by ELISA at 1:62500. HNRPA3 antibody can be used for detection of HNRPA3 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Heart at 4...



Antibody used in IHC on Human Liver at 4...



Antibody used in WB on Human Jurkat cell...