


Product Datasheet

SEPSECS Antibody (orb1241949)

Description	SEPSECS Antibody	
Species/Host	Rabbit	
Reactivity	Canine, Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	Antibody used in WB on Human Jurkat 0.2-...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLA/LP.	
Target	SEPSECS	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SLA/LP antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in IHC on Human kidney at ...
Note	For research use only	
Application notes	SLA/LP antibody can be used for detection of SLA/LP by ELISA at 1:312500. SLA/LP antibody can be used for detection of SLA/LP by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	49 kDa	
Uniprot ID	A1A4F3	
NCBI	NP_722547	
Dilution Range	SLA/LP antibody can be used for detection of SLA/LP by ELISA at 1:312500. SLA/LP antibody can be used for detection of SLA/LP by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	