

Product Datasheet

MATR3 Antibody (orb1242017)

Description	MATR3 Antibody	
Species/Host	Rabbit	
Reactivity	Human	Antibody used in IHC on Human Heart at 4...
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	Antibody used in IHC on Human Muscle at ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MATR3.	
Target	MATR3	Antibody used in WB on Human Jurkat cell...
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store MATR3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	MATR3 antibody can be used for detection of MATR3 by ELISA at 1:312500. MATR3 antibody can be used for detection of MATR3 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	93 kDa, 95 kDa	
Uniprot ID	P43243	
NCBI	NP_061322	
Dilution Range	MATR3 antibody can be used for detection of MATR3 by ELISA at 1:312500. MATR3 antibody can be used for detection of MATR3 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	