


Product Datasheet

RBM28 Antibody (orb1242021)


Description

RBM28 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RBM28.
Target	RBM28
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RBM28 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	RBM28 antibody can be used for detection of RBM28 by ELISA at 1:312500. RBM28 antibody can be used for detection of RBM28 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	83 kDa
Uniprot ID	Q9NW13
NCBI	NP_060547
Dilution Range	RBM28 antibody can be used for detection of RBM28 by ELISA at 1:312500. RBM28 antibody can be used for detection of RBM28 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Skin at 4....



Antibody used in IHC on Human Heart at 4...



Antibody used in WB on Human HepG2 at 0....