

Product Datasheet

RBM22 Antibody (orb1242022)

Description

RBM22 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RBM22.

Target

RBM22

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store RBM22 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

RBM22 antibody can be used for detection of RBM22 by ELISA at 1:1562500. RBM22 antibody can be used for detection of RBM22 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

46 kDa

Uniprot ID

[Q9NW64](#)

NCBI


[NP_060517](#)

Dilution Range


RBM22 antibody can be used for detection of RBM22 by ELISA at 1:1562500. RBM22 antibody can be used for detection of RBM22 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Heart at 4...



Antibody used in IHC on Human Muscle at ...



Antibody used in WB on Human Jurkat 0.2-...