
Product Datasheet

PRPF6 Antibody (orb1242054)

Description

PRPF6 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

 Antibody
used in WB
on
Transfected
293T ...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PRPF6.

Target

PRPF6

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PRPF6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PRPF6 antibody can be used for detection of PRPF6 by ELISA at 1:62500. PRPF6 antibody can be used for detection of PRPF6 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

104 kDa

Uniprot ID
[O94906](#)
NCBI
[NP_036601](#)
Dilution Range

PRPF6 antibody can be used for detection of PRPF6 by ELISA at 1:62500. PRPF6 antibody can be used for detection of PRPF6 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.