
Product Datasheet

PRPF6 Antibody (orb1242055)

Description	PRPF6 Antibody	
Species/Host	Rabbit	100
Reactivity	Human, Mouse, Rat	100
Conjugation	Unconjugated	100
Tested Applications	ELISA, WB	100
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PRPF6.	Antibody used in WB on Transfected 293T ...
Target	PRPF6	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PRPF6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	PRPF6 antibody can be used for detection of PRPF6 by ELISA at 1:62500. PRPF6 antibody can be used for detection of PRPF6 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	104 kDa, 47 kDa	
Uniprot ID	O94906	
NCBI	NP_036601	
Dilution Range	PRPF6 antibody can be used for detection of PRPF6 by ELISA at 1:62500. PRPF6 antibody can be used for detection of PRPF6 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	