
Product Datasheet

IGF2BP1 Antibody (orb1242080)

Description

IGF2BP1 Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human IGF2BP1.

Target

IGF2BP1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store IGF2BP1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

IGF2BP1 antibody can be used for detection of IGF2BP1 by ELISA at 1:312500. IGF2BP1 antibody can be used for detection of IGF2BP1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

63 kDa


Uniprot ID
[Q9NZI8](#)
NCBI
[NP_006537](#)
Dilution Range

IGF2BP1 antibody can be used for detection of IGF2BP1 by ELISA at 1:312500. IGF2BP1 antibody can be used for detection of IGF2BP1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

IGFBP1



 Brown: IGFBP1 Blue: Nuclei
See HC1 Data and Customer Feedback for more information

Antibody used in IHC on Human lung adeno...



Antibody used in WB on Human HepG2 at 0....

IGFBP1


 Brown: IGFBP1 Blue: Nuclei
See HC2 Data and Customer Feedback for more information

Antibody used in IHC on Human Lung at: 1...