
Product Datasheet

SF3B4 Antibody (orb1242106)

Description

SF3B4 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SF3B4.

Target

SF3B4

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SF3B4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SF3B4 antibody can be used for detection of SF3B4 by ELISA at 1:1562500. SF3B4 antibody can be used for detection of SF3B4 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

47 kDa

Uniprot ID
[Q15427](#)
NCBI
[NP_005841](#)
Dilution Range


SF3B4 antibody can be used for detection of SF3B4 by ELISA at 1:1562500. SF3B4 antibody can be used for detection of SF3B4 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human HepG2 at 1....



Antibody used in IHC on Human Intestine ...