
Product Datasheet

EIF2S1 Antibody (orb1242138)

Description

EIF2S1 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human EIF2S1.

Target

EIF2S1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store EIF2S1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

EIF2S1 antibody can be used for detection of EIF2S1 by ELISA at 1:2500. EIF2S1 antibody can be used for detection of EIF2S1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW


35 kDa

Uniprot ID
[P05198](#)
NCBI
[NP_004085](#)
Dilution Range


EIF2S1 antibody can be used for detection of EIF2S1 by ELISA at 1:2500. EIF2S1 antibody can be used for detection of EIF2S1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Lung at 4...



Antibody used in IHC on human Kidney at ...



Antibody used in WB on Human HepG2 at 1....