


Product Datasheet

RBM4 Antibody (orb1242166)

Description	RBM4 Antibody	RBM4
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	Antibody used in WB on Hum. Adult Placen...
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RBM4.	
Target	RBM4	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RBM4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Human MCF7 cells ...
Note	For research use only	
Application notes	RBM4 antibody can be used for detection of RBM4 by ELISA at 1:62500. RBM4 antibody can be used for detection of RBM4 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	40 kDa	
Uniprot ID	Q96LT9	
NCBI	NP_002887	
Dilution Range	RBM4 antibody can be used for detection of RBM4 by ELISA at 1:62500. RBM4 antibody can be used for detection of RBM4 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	