


Product Datasheet

NOVA2 Antibody (orb1242175)

Description	NOVA2 Antibody	NOVA2
Species/Host	Rabbit	
Reactivity	Human, Zebrafish	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	<small>See Immunoblot Data and Customer Feedback for more information</small>
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NOVA2.	Antibody used in WB on Human neural tube...
Target	NOVA2	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store NOVA2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Human HepG2 at 1....
Note	For research use only	
Application notes	NOVA2 antibody can be used for detection of NOVA2 by ELISA at 1:1562500. NOVA2 antibody can be used for detection of NOVA2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	54 kDa	
Uniprot ID	Q9UNW9	
NCBI	NP_002507	
Dilution Range	NOVA2 antibody can be used for detection of NOVA2 by ELISA at 1:1562500. NOVA2 antibody can be used for detection of NOVA2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	