


Product Datasheet

DAZ2 Antibody (orb1242229)

Description	DAZ2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human DAZ2.
Target	DAZ2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store DAZ2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	DAZ2 antibody can be used for detection of DAZ2 by ELISA at 1:62500. DAZ2 antibody can be used for detection of DAZ2 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	59 kDa, 65 kDa, 65 kDa, 41 kDa, 65 kDa, 65 kDa, 63 kDa, 50 kDa, 44 kDa
Uniprot ID	Q2KHN6
NCBI	NP_001005785
Dilution Range	DAZ2 antibody can be used for detection of DAZ2 by ELISA at 1:62500. DAZ2 antibody can be used for detection of DAZ2 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Skin at 4....



Antibody used in WB on Human Daudi cells...