

Product Datasheet

RPL32 Antibody (orb1242242)

Description	RPL32 Antibody	
Species/Host	Rabbit	
Reactivity	Canine, Human, Mouse, Rat	<p>Antibody used in WB on Human HepG2 at 2....</p>
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RPL32.	
Target	RPL32	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RPL32 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	RPL32 antibody can be used for detection of RPL32 by ELISA at 1:312500. RPL32 antibody can be used for detection of RPL32 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	15 kDa, 16 kDa, 16 kDa, 16 kDa, 16 kDa, 16 kDa, 16 kDa, 16 kDa	
Uniprot ID	P62910	
NCBI	NP_000985	
Dilution Range	RPL32 antibody can be used for detection of RPL32 by ELISA at 1:312500. RPL32 antibody can be used for detection of RPL32 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	