
Product Datasheet

ZBTB9 Antibody (orb1242303)

Description	ZBTB9 Antibody	kDa
Species/Host	Rabbit	40
Reactivity	Human	30
Conjugation	Unconjugated	20
Tested Applications	ELISA, WB	Antibody used in WB on Transfected 293T ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZBTB9.	
Target	ZBTB9	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ZBTB9 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	ZBTB9 antibody can be used for detection of ZBTB9 by ELISA at 1:312500. ZBTB9 antibody can be used for detection of ZBTB9 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	51 kDa	
Uniprot ID	Q96C00	
NCBI	NP_689948	
Dilution Range	ZBTB9 antibody can be used for detection of ZBTB9 by ELISA at 1:312500. ZBTB9 antibody can be used for detection of ZBTB9 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	