
Product Datasheet

KDM3B Antibody (orb1242420)

Description	KDM3B Antibody	
Species/Host	Rabbit	
Reactivity	Canine, Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human JMJD1B.	Antibody used in WB on Human Thymus at 2...
Target	KDM3B	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store JMJD1B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	JMJD1B antibody can be used for detection of JMJD1B by ELISA at 1:12500. JMJD1B antibody can be used for detection of JMJD1B by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	192 kDa	
Uniprot ID	Q7LBC6	
NCBI	NP_057688	
Dilution Range	JMJD1B antibody can be used for detection of JMJD1B by ELISA at 1:12500. JMJD1B antibody can be used for detection of JMJD1B by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	