

Product Datasheet

SAP30BP Antibody (orb1242458)

Description

SAP30BP Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

 Antibody
used in WB
on Human
Muscle at 1...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SAP30BP.

Target

SAP30BP

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SAP30BP antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SAP30BP antibody can be used for detection of SAP30BP by ELISA at 1:1562500. SAP30BP antibody can be used for detection of SAP30BP by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

34 kDa

Uniprot ID
[Q9UHR5](#)
NCBI
[NP_037392](#)
Dilution Range

SAP30BP antibody can be used for detection of SAP30BP by ELISA at 1:1562500. SAP30BP antibody can be used for detection of SAP30BP by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.