
Product Datasheet

PSMC5 Antibody (orb1242631)

Description

PSMC5 Antibody

Species/Host

Rabbit

Reactivity

C. elegans, Canine, Drosophila, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested Applications

ELISA, WB

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PSMC5.

Target

PSMC5

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PSMC5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PSMC5 antibody can be used for detection of PSMC5 by ELISA at 1:312500. PSMC5 antibody can be used for detection of PSMC5 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

46 kDa

Uniprot ID
[P62195](#)
NCBI
[NP_002796](#)
Dilution Range

PSMC5 antibody can be used for detection of PSMC5 by ELISA at 1:312500. PSMC5 antibody can be used for detection of PSMC5 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

PSMC5



Antibody used in WB on Mouse Brain at 2 ...

PSMC5

