

Product Datasheet

FOXO3 Antibody (orb1242669)


Description	FOXO3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human FOXO3.
Target	FOXO3
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store FOXO3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	FOXO3 antibody can be used for detection of FOXO3 by ELISA at 1:312500. FOXO3 antibody can be used for detection of FOXO3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	71 kDa
Uniprot ID	O43524
NCBI	NP_001446
Dilution Range	FOXO3 antibody can be used for detection of FOXO3 by ELISA at 1:312500. FOXO3 antibody can be used for detection of FOXO3 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Transfected 293T ...



Antibody used in WB on Hum. Fetal Heart ...



Antibody used in WB on Hum. Fetal Lung a...