
Product Datasheet

Tbpl2 Antibody (orb1242738)

Description	Tbpl2 Antibody	SP2/0
Species/Host	Rabbit	SP2/0
Reactivity	Mouse	SP2/0
Conjugation	Unconjugated	SP2/0
Tested Applications	ELISA, WB	Antibody used in WB on Mouse SP2/0 at 5...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse TRF3.	
Target	Tbpl2	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store TRF3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	TRF3 antibody can be used for detection of TRF3 by ELISA at 1:1562500. TRF3 antibody can be used for detection of TRF3 by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	39 kDa	
Uniprot ID	Q6SJ95	
NCBI	NP_951014	
Dilution Range	TRF3 antibody can be used for detection of TRF3 by ELISA at 1:1562500. TRF3 antibody can be used for detection of TRF3 by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	