

Product Datasheet

Mri1 Antibody (orb1242785)

Description	Mri1 Antibody	
Species/Host	Rabbit	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Mouse NIH-3T3 at ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human 2410018C20RIK.	
Target	Mri1	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store 2410018C20RIK antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	2410018C20RIK antibody can be used for detection of 2410018C20RIK by ELISA at 1:1562500. 2410018C20RIK antibody can be used for detection of 2410018C20RIK by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	41 kDa	
Uniprot ID	Q9CQT1	
NCBI	NP_080699	
Dilution Range	2410018C20RIK antibody can be used for detection of 2410018C20RIK by ELISA at 1:1562500. 2410018C20RIK antibody can be used for detection of 2410018C20RIK by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	