

Product Datasheet

ASCL2 Antibody (orb1242856)

Description	ASCL2 Antibody	
Species/Host	Rabbit	
Reactivity	Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human Small intes...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse ASCL2.	
Target	ASCL2	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ASCL2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	ASCL2 antibody can be used for detection of ASCL2 by ELISA at 1:312500. ASCL2 antibody can be used for detection of ASCL2 by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	29 kDa	
Uniprot ID	O35885	
NCBI	NP_032580	
Dilution Range	ASCL2 antibody can be used for detection of ASCL2 by ELISA at 1:312500. ASCL2 antibody can be used for detection of ASCL2 by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	