
Product Datasheet

SP140L Antibody (orb1243002)

Description	SP140L Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human Jurkat 2.5 ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human LOC93349.	
Target	SP140L	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store LOC93349 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	LOC93349 antibody can be used for detection of LOC93349 by ELISA at 1:1562500. LOC93349 antibody can be used for detection of LOC93349 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	39 kDa	
Uniprot ID	Q9H930	
NCBI	NP_612411	
Dilution Range	LOC93349 antibody can be used for detection of LOC93349 by ELISA at 1:1562500. LOC93349 antibody can be used for detection of LOC93349 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	