

Product Datasheet

ZNF416 Antibody (orb1243053)

Description	ZNF416 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human 293T at 1.2...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZNF416.	
Target	ZNF416	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ZNF416 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	ZNF416 antibody can be used for detection of ZNF416 by ELISA at 1:1562500. ZNF416 antibody can be used for detection of ZNF416 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	67 kDa	
Uniprot ID	Q9BWM5	
NCBI	NP_060349	
Dilution Range	ZNF416 antibody can be used for detection of ZNF416 by ELISA at 1:1562500. ZNF416 antibody can be used for detection of ZNF416 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	