
Product Datasheet

PIAS2 Antibody (orb1243100)

Description	PIAS2 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	91kDa
Conjugation	Unconjugated	47kDa
Tested Applications	ELISA, WB	32kDa
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PIAS2.	Antibody used in WB on Human HepG2 at 0....
Target	PIAS2	PIAS2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	91kDa
Form/Appearance	Liquid	65kDa
Concentration	batch dependent	40kDa
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PIAS2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	29kDa
Note	For research use only	22kDa
Application notes	PIAS2 antibody can be used for detection of PIAS2 by ELISA at 1:312500. PIAS2 antibody can be used for detection of PIAS2 by western blot at 0.6 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	Antibody used in WB on Human Hela at 1 u...
Clonality	Polyclonal	PIAS2
MW	68 kDa	91kDa
Uniprot ID	O75928	65kDa
NCBI	NP_004662	40kDa
Dilution Range	PIAS2 antibody can be used for detection of PIAS2 by ELISA at 1:312500. PIAS2 antibody can be used for detection of PIAS2 by western blot at 0.6 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	29kDa
Expiration Date	12 months from date of receipt.	22kDa