
Product Datasheet

ARIH2 Antibody (orb1243351)

Description

ARIH2 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ARIH2.

Target

ARIH2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ARIH2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ARIH2 antibody can be used for detection of ARIH2 by ELISA at 1:62500. ARIH2 antibody can be used for detection of ARIH2 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW


58 kDa

Uniprot ID
[O95376](#)
NCBI
[NP_006312](#)
Dilution Range

ARIH2 antibody can be used for detection of ARIH2 by ELISA at 1:62500. ARIH2 antibody can be used for detection of ARIH2 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Lung.



Antibody used in WB on Human Jurkat 2.5 ...



Antibody used in WB on Jurkat 1.25 µg/ml...