

Product Datasheet

TRIM13 Antibody (orb1243421)

Description

TRIM13 Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RFP2.

Target

TRIM13

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store RFP2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

RFP2 antibody can be used for detection of RFP2 by ELISA at 1:62500. RFP2 antibody can be used for detection of RFP2 by western blot at 0.125 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

47 kDa, 47 kDa, 47 kDa, 47 kDa

Uniprot ID

[Q5UBW0](#)

NCBI


[NP_001007279](#)

Dilution Range


RFP2 antibody can be used for detection of RFP2 by ELISA at 1:62500. RFP2 antibody can be used for detection of RFP2 by western blot at 0.125 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in IHC on Human Lung.



Antibody used in WB on Human Lung at 0.1...