
Product Datasheet

KMT2B Antibody (orb1243505)

Description

KMT2B Antibody

Species/Host

Rabbit

Reactivity


Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications


Antibody used in WB on Human HepG2 cells...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MLL4.

Target

KMT2B

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store MLL4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

MLL4 antibody can be used for detection of MLL4 by ELISA at 1:12500. MLL4 antibody can be used for detection of MLL4 by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

64 kDa, 294 kDa

Uniprot ID
[Q9UMN6](#)
NCBI
[Q9UMN6](#)
Dilution Range

MLL4 antibody can be used for detection of MLL4 by ELISA at 1:12500. MLL4 antibody can be used for detection of MLL4 by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.