

Product Datasheet

SUPT5H Antibody (orb1243583)

Description

SUPT5H Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SUPT5H.

Target

SUPT5H

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SUPT5H antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SUPT5H antibody can be used for detection of SUPT5H by ELISA at 1:62500. SUPT5H antibody can be used for detection of SUPT5H by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

121 kDa

Uniprot ID

[O00267](#)

NCBI


[NP_003160](#)

Dilution Range


SUPT5H antibody can be used for detection of SUPT5H by ELISA at 1:62500. SUPT5H antibody can be used for detection of SUPT5H by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Pancreas.



Antibody used in IHC on Human kidney.



Antibody used in WB on Human HepG2 at 0....