
Product Datasheet

SOX18 Antibody (orb1243651)

Description

SOX18 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SOX18.

Target

SOX18

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SOX18 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SOX18 antibody can be used for detection of SOX18 by ELISA at 1:62500. SOX18 antibody can be used for detection of SOX18 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW


41 kDa

Uniprot ID
[P35713](#)
NCBI
[NP_060889](#)
Dilution Range

SOX18 antibody can be used for detection of SOX18 by ELISA at 1:62500. SOX18 antibody can be used for detection of SOX18 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human kidney at ...



Antibody used in WB on Human heart at 0....



Antibody used in IF on Human lymphatic e...