
Product Datasheet

GTF2A1L Antibody (orb1243703)

Description

GTF2A1L Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, WB

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ALF.

Target

GTF2A1L

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ALF antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ALF antibody can be used for detection of ALF by ELISA at 1:1562500. ALF antibody can be used for detection of ALF by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

50 kDa

Uniprot ID
[Q9UNN4](#)
NCBI
[NP_751946](#)
Dilution Range

ALF antibody can be used for detection of ALF by ELISA at 1:1562500. ALF antibody can be used for detection of ALF by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

Antibody used in WB on Human Jurkat 0.2-...