
Product Datasheet

POLR2B Antibody (orb1243819)

Description

POLR2B Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody used in WB on Human HepG2 at 5

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human POLR2B.

Target

POLR2B

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store POLR2B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

POLR2B antibody can be used for detection of POLR2B by ELISA at 1:312500. POLR2B antibody can be used for detection of POLR2B by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

134 kDa

Uniprot ID
[P30876](#)
NCBI
[NP_000929](#)
Dilution Range

POLR2B antibody can be used for detection of POLR2B by ELISA at 1:312500. POLR2B antibody can be used for detection of POLR2B by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.