

Product Datasheet

POLR2B Antibody (orb1243820)

Description	POLR2B Antibody	
Species/Host	Rabbit	
Reactivity	Canine, Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human Muscle at 0...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human POLR2B.	
Target	POLR2B	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store POLR2B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	POLR2B antibody can be used for detection of POLR2B by ELISA at 1:1562500. POLR2B antibody can be used for detection of POLR2B by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	134 kDa	
Uniprot ID	P30876	
NCBI	NP_000929	
Dilution Range	POLR2B antibody can be used for detection of POLR2B by ELISA at 1:1562500. POLR2B antibody can be used for detection of POLR2B by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	