

Product Datasheet

OTX2 Antibody (orb1243827)

Description

OTX2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human OTX2.

Target

OTX2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store OTX2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

OTX2 antibody can be used for detection of OTX2 by ELISA at 1:312500. OTX2 antibody can be used for detection of OTX2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

32 kDa

Uniprot ID

[P32243](#)

NCBI


[NP_068374](#)

Dilution Range


OTX2 antibody can be used for detection of OTX2 by ELISA at 1:312500. OTX2 antibody can be used for detection of OTX2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Stomach.



Antibody used in IHC on Human Heart.



Antibody used in WB on Human Jurkat 0.2-...