

Product Datasheet

TAB2 Antibody (orb1243836)

Description

TAB2 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MAP3K7IP2.

Target

TAB2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store MAP3K7IP2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

MAP3K7IP2 antibody can be used for detection of MAP3K7IP2 by ELISA at 1:312500. MAP3K7IP2 antibody can be used for detection of MAP3K7IP2 by western blot at 2.0/mL.

Clonality

Polyclonal

MW

76 kDa

Uniprot ID

[Q9NYJ8](#)

NCBI


[NP_055908](#)

Dilution Range


MAP3K7IP2 antibody can be used for detection of MAP3K7IP2 by ELISA at 1:312500. MAP3K7IP2 antibody can be used for detection of MAP3K7IP2 by western blot at 2.0/mL.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Stomach.



Antibody used in IHC on Human Spermatozoa.



Antibody used in WB on Human HepG2 cells.