


Product Datasheet

ZFP1 Antibody (orb1243902)

Description	ZFP1 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZFP1.	Antibody used in WB on Human HepG2 at 2....
Target	ZFP1	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ZFP1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	ZFP1 antibody can be used for detection of ZFP1 by ELISA at 1:1562500. ZFP1 antibody can be used for detection of ZFP1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	41 kDa	
Uniprot ID	Q6P2D0	
NCBI	NP_710155	
Dilution Range	ZFP1 antibody can be used for detection of ZFP1 by ELISA at 1:1562500. ZFP1 antibody can be used for detection of ZFP1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	