

Product Datasheet

EBF1 Antibody (orb1244037)

Description	EBF1 Antibody	
Species/Host	Rabbit	
Reactivity	C. elegans, Drosophila, Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human brain at 0...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human EBF1.	
Target	EBF1	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store EBF1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	EBF1 antibody can be used for detection of EBF1 by ELISA at 1:312500. EBF1 antibody can be used for detection of EBF1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	64 kDa	
Uniprot ID	Q9UH73	
NCBI	NP_076870	
Dilution Range	EBF1 antibody can be used for detection of EBF1 by ELISA at 1:312500. EBF1 antibody can be used for detection of EBF1 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	