
Product Datasheet

OLIG2 Antibody (orb1244060)

Description

OLIG2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, WB

Applications


Antibody used in WB on human LCL and mou...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human OLIG2.

Target

OLIG2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store OLIG2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

OLIG2 antibody can be used for detection of OLIG2 by ELISA at 1:62500. OLIG2 antibody can be used for detection of OLIG2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.



Antibody used in WB on Human 293T at 0.2...

Clonality

Polyclonal

MW

32 kDa

Uniprot ID
[Q13516](https://www.uniprot.org/uniprot/Q13516)
NCBI
[NP_005797](https://www.ncbi.nlm.nih.gov/nuccore/NP_005797)
Dilution Range

OLIG2 antibody can be used for detection of OLIG2 by ELISA at 1:62500. OLIG2 antibody can be used for detection of OLIG2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.



Antibody used in IHC on human LCL and mo...

Expiration Date

12 months from date of receipt.