
Product Datasheet

E2F4 Antibody (orb1244097)

Description

E2F4 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human E2F4.

Target

E2F4

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store E2F4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

E2F4 antibody can be used for detection of E2F4 by ELISA at 1:312500. E2F4 antibody can be used for detection of E2F4 by western blot at 2.0/mL.

Clonality

Polyclonal

MW


44 kDa

Uniprot ID
[Q16254](#)
NCBI
[NP_001941](#)
Dilution Range


E2F4 antibody can be used for detection of E2F4 by ELISA at 1:312500. E2F4 antibody can be used for detection of E2F4 by western blot at 2.0/mL.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Stomach.



Antibody used in IHC on Human Intestine.



Antibody used in WB on Human HepG2 at 0....