
Product Datasheet

NFIC Antibody (orb1244168)

| Description | NFIC Antibody | 800_ |
|----------------------------|--|---|
| Species/Host | Rabbit | 400_ |
| Reactivity | Human | 300_ |
| Conjugation | Unconjugated | 200_ |
| Tested Applications | ELISA, WB | Antibody used in WB on Human 721_B at 0.... |
| Immunogen | Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NFIC. | |
| Target | NFIC | |
| Preservatives | Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. | |
| Form/Appearance | Liquid | |
| Concentration | batch dependent | |
| Storage | For short periods of storage (days) store at 4°C. For longer periods of storage, store NFIC antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles. | |
| Note | For research use only | |
| Application notes | NFIC antibody can be used for detection of NFIC by ELISA at 1:62500. NFIC antibody can be used for detection of NFIC by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000. | |
| Clonality | Polyclonal | |
| MW | 55 kDa | |
| Uniprot ID | P08651 | |
| NCBI | NP_995315 | |
| Dilution Range | NFIC antibody can be used for detection of NFIC by ELISA at 1:62500. NFIC antibody can be used for detection of NFIC by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000. | |
| Expiration Date | 12 months from date of receipt. | |